

Annexin A4

Catalog No: #AP80760



Package Size: #AP80760-1 50ug #AP80760-2 100ug #AP80760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	35-beta calcimedlin,Annexin IV,Annexin-4,Carbohydrate-binding protein p33/p41,Chromobindin-4,Endonexin I,Lipocortin IV,P32.5,PP4-X,Placental anticoagulant protein II,PAP-II,Protein II,ANX4
Accession No.	Uniprot ID:P13214
Uniprot	P13214
GeneID	281625;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ANXA4

Note: This product is for in vitro research use only