

Ubiquitin-like protein ATG12

Catalog No: #AP80764



Package Size: #AP80764-1 50ug #AP80764-2 100ug #AP80764-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-like protein ATG12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Autophagy-related protein 12,APG12-like,Apg12,Apg12l
Accession No.	Uniprot ID:Q9CQY1
Uniprot	Q9CQY1
GeneID	67526;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atg12

Note: This product is for in vitro research use only