

DNA-(apurinic or apyrimidinic site) lyase

Catalog No: #AP80766



Package Size: #AP80766-1 50ug #AP80766-2 100ug #AP80766-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA-(apurinic or apyrimidinic site) lyase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	APEX nuclease, APEN, Apurinic-apyrimidinic endonuclease 1, AP endonuclease 1, REF-1, Redox factor-1, APE, APEX, BAP1, REF1
Accession No.	Uniprot ID:P23196
Uniprot	P23196
GeneID	281630;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APEX1

Note: This product is for in vitro research use only