

Apolipoprotein C-IV

Catalog No: #AP80768



Package Size: #AP80768-1 50ug #AP80768-2 100ug #AP80768-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-IV
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apo-CIV, Apolipoprotein C2-linked, ACL, Apolipoprotein C4, Acl
Accession No.	Uniprot ID:Q61268
Uniprot	Q61268
GeneID	11425;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Apoc4

Note: This product is for in vitro research use only