

a-actinin Monoclonal Antibody(4A12)

Catalog No: #44139

Package Size: #44139-1 50ul #44139-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	a-actinin Monoclonal Antibody(4A12)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A12
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	a-actinin Monoclonal Antibody(4A12) detects endogenous levels of ACTN1
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	a-actinin
Other Names	ACTN1
Accession No.	Swiss-Prot:P12814Gene ID:87
Uniprot	P12814
GeneID	87;
Calculated MW	103kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using a-actinin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using α -actinin Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only