

ATP synthase-coupling factor 6, mitochondrial

Catalog No: #AP80773

Package Size: #AP80773-1 50ug #AP80773-2 100ug #AP80773-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP synthase-coupling factor 6, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ATPase subunit F6
Accession No.	Uniprot ID:P02721
Uniprot	P02721
GeneID	282690;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATP5J

Note: This product is for in vitro research use only