

Axin-1

Catalog No: #AP80775



Package Size: #AP80775-1 50ug #AP80775-2 100ug #AP80775-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Axin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Axis inhibition protein 1,rAxin,Axin
Accession No.	Uniprot ID:O70239
Uniprot	O70239
GeneID	79257;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Axin1

Note: This product is for in vitro research use only