

Sodium/potassium-transporting ATPase subunit alpha-3

Catalog No: #AP80777

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80777-1 50ug #AP80777-2 100ug #AP80777-3 1mg

Description

Product Name	Sodium/potassium-transporting ATPase subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Na(+)/K(+) ATPase alpha-3 subunit,Na(+)/K(+) ATPase alpha(III) subunit,Sodium pump subunit alpha-3
Accession No.	Uniprot ID:P24798
Uniprot	P24798
GeneID	396467;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP1A3

Note: This product is for in vitro research use only