

Asialoglycoprotein receptor 1

Catalog No: #AP80781



Package Size: #AP80781-1 50ug #AP80781-2 100ug #AP80781-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Asialoglycoprotein receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ASGP-R 1,C-type lectin domain family 4 member H1,Hepatic lectin H1,HL-1,CLEC4H1
Accession No.	Uniprot ID:P07306
Uniprot	P07306
GeneID	432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ASGR1

Note: This product is for in vitro research use only