

Ataxin-10

Catalog No: #AP80782



Package Size: #AP80782-1 50ug #AP80782-2 100ug #AP80782-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ataxin-10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neuronal beta-catenin-like protein, Spinocerebellar ataxia type 10 protein homolog, Sca10
Accession No.	Uniprot ID:Q9ER24
Uniprot	Q9ER24
GeneID	170821;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atxn10

Note: This product is for in vitro research use only