

Cyclic AMP-dependent transcription factor ATF-4

Catalog No: #AP80783



Package Size: #AP80783-1 50ug #AP80783-2 100ug #AP80783-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	cAMP-dependent transcription factor ATF-4,Activating transcription factor 4,rATF-4
Accession No.	Uniprot ID:Q9ES19
Uniprot	Q9ES19
GeneID	79255;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atf4

Note: This product is for in vitro research use only