

ATP synthase-coupling factor 6, mitochondrial

Catalog No: #AP80784



Package Size: #AP80784-1 50ug #AP80784-2 100ug #AP80784-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP synthase-coupling factor 6, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATPase subunit F6,ATP5A,ATPM
Accession No.	Uniprot ID:P18859
Uniprot	P18859
GeneID	522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP5J

Note: This product is for in vitro research use only