

Gamma-secretase subunit APH-1A

Catalog No: #AP80786



Package Size: #AP80786-1 50ug #AP80786-2 100ug #AP80786-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Gamma-secretase subunit APH-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APH-1a,Aph-1alpha,Presenilin-stabilization factor,PSF,CGI-78,UNQ579/PRO1141
Accession No.	Uniprot ID:Q96BI3
Uniprot	Q96BI3
GeneID	51107;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APH1A

Note: This product is for in vitro research use only