

Aryl hydrocarbon receptor nuclear translocator

Catalog No: #AP80788



Package Size: #AP80788-1 50ug #AP80788-2 100ug #AP80788-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aryl hydrocarbon receptor nuclear translocator
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ARNT protein,Dioxin receptor, nuclear translocator,Hypoxia-inducible factor 1-beta,HIF-1-beta
Accession No.	Uniprot ID:P41739
Uniprot	P41739
GeneID	25242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Arnt

Note: This product is for in vitro research use only