

## Plasma membrane calcium-transporting ATPase 2

Catalog No: #AP80789



Package Size: #AP80789-1 50ug #AP80789-2 100ug #AP80789-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Plasma membrane calcium-transporting ATPase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, PMCA2
Accession No.	Uniprot ID:Q01814
Uniprot	Q01814
GeneID	491;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ATP2B2

Note: This product is for in vitro research use only