

## IL-8 Mouse Monoclonal Antibody(8B1)

Catalog No: #44141

Package Size: #44141-1 50ul #44141-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IL-8 Mouse Monoclonal Antibody(8B1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8B1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	IL-8 Mouse Monoclonal Antibody(8B1) detects endogenous levels of IL8
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	IL-8
Other Names	IL8
Accession No.	Swiss-Prot:P10145Gene ID:3576
Uniprot	P10145
GeneID	3576;
Calculated MW	11kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using IL-8 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using IL-8 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only