

Cyclic AMP-dependent transcription factor ATF-4

Catalog No: #AP80790



Package Size: #AP80790-1 50ug #AP80790-2 100ug #AP80790-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	cAMP-dependent transcription factor ATF-4,Activating transcription factor 4,Cyclic AMP-responsive element-binding protein 2,CREB-2,DNA-binding protein TAXREB67,Tax-responsive enhancer element-binding protein 67,TaxREB67,CREB2,TXREB
Accession No.	Uniprot ID:P18848
Uniprot	P18848
GeneID	468;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATF4

Note: This product is for in vitro research use only