

AP-1 complex subunit mu-1

Catalog No: #AP80795

Package Size: #AP80795-1 50ug #AP80795-2 100ug #AP80795-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | AP-1 complex subunit mu-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | AP-mu chain family member mu1A,Adaptor protein complex AP-1 subunit mu-1,Adaptor-related protein complex 1 subunit mu-1,Clathrin assembly protein complex 1 mu-1 medium chain 1,Clathrin coat assembly protein AP47,Clathrin coat-associated protein AP47,Golgi adaptor HA1/AP1 adaptin mu-1 subunit,Mu-adaptin 1,Mu1A-adaptin,Cltnm |
| Accession No. | Uniprot ID:P35585 |
| Uniprot | P35585 |
| GeneID | 11767; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ap1m1

Note: This product is for in vitro research use only