

## AP-2 complex subunit beta

Catalog No: #AP80796



Package Size: #AP80796-1 50ug #AP80796-2 100ug #AP80796-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AP-2 complex subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AP105B,Adaptor protein complex AP-2 subunit beta,Adaptor-related protein complex 2 subunit beta,Beta-2-adaptin,Beta-adaptin,Clathrin assembly protein complex 2 beta large chain,Plasma membrane adaptor HA2/AP2 adaptin beta subunit,Clapp1
Accession No.	Uniprot ID:Q9DBG3
Uniprot	Q9DBG3
GeneID	71770;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ap2b1

Note: This product is for in vitro research use only