

Annexin A2

Catalog No: #AP80800



Package Size: #AP80800-1 50ug #AP80800-2 100ug #AP80800-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Annexin A2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Annexin II,Annexin-2,Calpactin I heavy chain,Calpactin-1 heavy chain,Chromobindin-8,Lipocortin II,Placental anticoagulant protein IV,PAP-IV,Protein I,p36,ANX2 |
| Accession No. | Uniprot ID:P04272 |
| Uniprot | P04272 |
| GeneID | 282689; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ANXA2

Note: This product is for in vitro research use only