

## Sarcoplasmic/endoplasmic reticulum calcium ATPase

1

Catalog No: #AP80806

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80806-1 50ug #AP80806-2 100ug #AP80806-3 1mg

## Description

Product Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Uniprot ID:Q8R429
Uniprot	Q8R429
GeneID	11937;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Atp2a1

Note: This product is for in vitro research use only