

# ATP synthase-coupling factor 6, mitochondrial

Catalog No: #AP80808



Package Size: #AP80808-1 50ug #AP80808-2 100ug #AP80808-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | ATP synthase-coupling factor 6, mitochondrial          |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms                                                     |
| Other Names        | ATPase subunit F6                                      |
| Accession No.      | Uniprot ID:P97450                                      |
| Uniprot            | P97450                                                 |
| GeneID             | 11957;                                                 |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Atp5j

Note: This product is for in vitro research use only