

HDAC1 Mouse Monoclonal Antibody(4E1)

Catalog No: #44143

Package Size: #44143-1 50ul #44143-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HDAC1 Mouse Monoclonal Antibody(4E1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4.00E+01
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	HDAC1 Mouse Monoclonal Antibody(4E1) detects endogenous levels of HDAC1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	HDAC1
Other Names	HDAC1
Accession No.	Swiss-Prot:Q13547Gene ID:3065
Uniprot	Q13547
GeneID	3065;
Calculated MW	62kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HDAC1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using HADC1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only