

Artemin

Catalog No: #AP80811



Package Size: #AP80811-1 50ug #AP80811-2 100ug #AP80811-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Artemin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Enovin,Neublastin,EVN
Accession No.	Uniprot ID:Q5T4W7
Uniprot	Q5T4W7
GeneID	9048;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ARTN

Note: This product is for in vitro research use only