

Sarcoplasmic/endoplasmic reticulum calcium ATPase

3



Catalog No: #AP80813

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80813-1 50ug #AP80813-2 100ug #AP80813-3 1mg

Description

Product Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SERCA3, Calcium pump 3
Accession No.	Uniprot ID:Q64518
Uniprot	Q64518
GeneID	53313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atp2a3

Note: This product is for in vitro research use only