

Plasma membrane calcium-transporting ATPase 2

Catalog No: #AP80815



Package Size: #AP80815-1 50ug #AP80815-2 100ug #AP80815-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Plasma membrane calcium-transporting ATPase 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2 |
| Accession No. | Uniprot ID:P11506 |
| Uniprot | P11506 |
| GeneID | 24215; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atp2b2

Note: This product is for in vitro research use only