

Aprataxin

Catalog No: #AP80816



Package Size: #AP80816-1 50ug #AP80816-2 100ug #AP80816-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aprataxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Forkhead-associated domain histidine triad-like protein,FHA-HIT
Accession No.	Uniprot ID:Q7YRZ2
Uniprot	Q7YRZ2
GeneID	359714;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APTX

Note: This product is for in vitro research use only