

Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1

Catalog No: #AP80821

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80821-1 50ug #AP80821-2 100ug #AP80821-3 1mg

Description

Product Name	Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B3Gal-T8,Core 1 O-glycan T-synthase,Core 1 UDP-galactose:N-acetylgalactosamine-alpha-R beta 1,3-galactosyltransferase 1,Beta-1,3-galactosyltransferase,Core 1 beta1,3-galactosyltransferase 1,C1GalT1
Accession No.	Uniprot ID:Q9NS00
Uniprot	Q9NS00
GeneID	56913;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C1GALT1

Note: This product is for in vitro research use only