

CD209 antigen-like protein A

Catalog No: #AP80822



Package Size: #AP80822-1 50ug #AP80822-2 100ug #AP80822-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD209 antigen-like protein A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dendritic cell-specific ICAM-3-grabbing non-integrin,DC-SIGN,CD209,Cire
Accession No.	Uniprot ID:Q91ZX1
Uniprot	Q91ZX1
GeneID	170786;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd209a

Note: This product is for in vitro research use only