

## Isoaspartyl peptidase/L-asparaginase

Catalog No: #AP80829



Package Size: #AP80829-1 50ug #AP80829-2 100ug #AP80829-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Isoaspartyl peptidase/L-asparaginase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Asparaginase-like protein 1, Beta-aspartyl-peptidase, Isoaspartyl dipeptidase, L-asparagine amidohydrolase
Accession No.	Uniprot ID:Q32LE5
Uniprot	Q32LE5
GeneID	767970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ASRGL1

Note: This product is for in vitro research use only