

Acetylserotonin O-methyltransferase

Catalog No: #AP80830



Package Size: #AP80830-1 50ug #AP80830-2 100ug #AP80830-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetylserotonin O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hydroxyindole O-methyltransferase,HIOMT
Accession No.	Uniprot ID:P10950
Uniprot	P10950
GeneID	281013;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ASMT

Note: This product is for in vitro research use only