

AT-rich interactive domain-containing protein 4B

Catalog No: #AP80832



Package Size: #AP80832-1 50ug #AP80832-2 100ug #AP80832-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AT-rich interactive domain-containing protein 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ARID domain-containing protein 4B, 180 kDa Sin3-associated polypeptide, Sin3-associated polypeptide p180, Breast cancer-associated antigen BRCAA1, Histone deacetylase complex subunit SAP180, Retinoblastoma-binding protein 1-like 1, BRCAA1, RBBP1L1, RBP1L1, SAP180
Accession No.	Uniprot ID:Q4LE39
Uniprot	Q4LE39
GeneID	51742;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ARID4B

Note: This product is for in vitro research use only