

## Beta-2-glycoprotein 1

Catalog No: #AP80834



Package Size: #AP80834-1 50ug #AP80834-2 100ug #AP80834-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-2-glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	APC inhibitor,Activated protein C-binding protein,Apolipoprotein H,Apo-H,Beta-2-glycoprotein 1,B2GPI,B2gp1
Accession No.	Uniprot ID:Q01339
Uniprot	Q01339
GeneID	11818;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ApoH

Note: This product is for in vitro research use only