

Aquaporin-2

Catalog No: #AP80838



Package Size: #AP80838-1 50ug #AP80838-2 100ug #AP80838-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aquaporin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AQP-2,ADH water channel,Aquaporin-CD,AQP-CD,Collecting duct water channel protein,WCH-CD,Water channel protein for renal collecting duct
Accession No.	Uniprot ID:P34080
Uniprot	P34080
GeneID	25386;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aqp2

Note: This product is for in vitro research use only