

Copper-transporting ATPase 1

Catalog No: #AP80840



Package Size: #AP80840-1 50ug #AP80840-2 100ug #AP80840-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Copper-transporting ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Copper pump 1, Menkes disease-associated protein, MC1, MNK
Accession No.	Uniprot ID:Q04656
Uniprot	Q04656
GeneID	538;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATP7A

Note: This product is for in vitro research use only