

Aminopeptidase N

Catalog No: #AP80846



Package Size: #AP80846-1 50ug #AP80846-2 100ug #AP80846-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aminopeptidase N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	AP-N,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Microsomal aminopeptidase,gp130,CD13
Accession No.	Uniprot ID:P15145
Uniprot	P15145
GeneID	397520;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ANPEP

Note: This product is for in vitro research use only