

Large proline-rich protein BAG6

Catalog No: #AP80849



Package Size: #AP80849-1 50ug #AP80849-2 100ug #AP80849-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Large proline-rich protein BAG6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BAG family molecular chaperone regulator 6,BCL2-associated athanogene 6,BAG-6,HLA-B-associated transcript 3,Protein Scythe,Bat3
Accession No.	Uniprot ID:Q9Z1R2
Uniprot	Q9Z1R2
GeneID	224727;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bag6

Note: This product is for in vitro research use only