

Autophagy-related protein 9A

Catalog No: #AP80852

Package Size: #AP80852-1 50ug #AP80852-2 100ug #AP80852-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Autophagy-related protein 9A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	APG9-like 1, Apg9l1
Accession No.	Uniprot ID:Q68FE2
Uniprot	Q68FE2
GeneID	245860;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atg9a

Note: This product is for in vitro research use only