

## DNA-(apurinic or apyrimidinic site) lyase

Catalog No: #AP80853



Package Size: #AP80853-1 50ug #AP80853-2 100ug #AP80853-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA-(apurinic or apyrimidinic site) lyase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	APEX nuclease, APEN, Apurinic-apyrimidinic endonuclease 1, AP endonuclease 1, REF-1, Redox factor-1, Ape, Apex, Ref1
Accession No.	Uniprot ID:P28352
Uniprot	P28352
GeneID	11792;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Apex1

Note: This product is for in vitro research use only