

Aprataxin and PNK-like factor

Catalog No: #AP80854



Package Size: #AP80854-1 50ug #AP80854-2 100ug #AP80854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aprataxin and PNK-like factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apurinic-apyrimidinic endonuclease APLF
Accession No.	Uniprot ID:Q9D842
Uniprot	Q9D842
GeneID	72103;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Aplf

Note: This product is for in vitro research use only