

Plasma membrane calcium-transporting ATPase 1

Catalog No: #AP80858



Package Size: #AP80858-1 50ug #AP80858-2 100ug #AP80858-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasma membrane calcium-transporting ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	PMCA1, Plasma membrane calcium ATPase isoform 1, Plasma membrane calcium pump isoform 1, PMCA1
Accession No.	Uniprot ID:Q98SH2
Uniprot	Q98SH2
GeneID	374244;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATP2B1

Note: This product is for in vitro research use only