

Cyclic AMP-dependent transcription factor ATF-4

Catalog No: #AP80859



Package Size: #AP80859-1 50ug #AP80859-2 100ug #AP80859-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	cAMP-dependent transcription factor ATF-4, Activating transcription factor 4
Accession No.	Uniprot ID:Q3ZCH6
Uniprot	Q3ZCH6
GeneID	509107;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATF4

Note: This product is for in vitro research use only