

Beta-galactosidase

Catalog No: #AP80861



Package Size: #AP80861-1 50ug #AP80861-2 100ug #AP80861-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-galactosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acid beta-galactosidase,Lactase,Elastin receptor 1,ELNR1
Accession No.	Uniprot ID:P16278
Uniprot	P16278
GeneID	2720;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLB1

Note: This product is for in vitro research use only