

Serine-protein kinase ATM

Catalog No: #AP80867



Package Size: #AP80867-1 50ug #AP80867-2 100ug #AP80867-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine-protein kinase ATM
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ataxia telangiectasia mutated homolog,A-T mutated homolog
Accession No.	Uniprot ID:Q62388
Uniprot	Q62388
GeneID	11920;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atm

Note: This product is for in vitro research use only