

Bone morphogenetic protein 5

Catalog No: #AP80870



Package Size: #AP80870-1 50ug #AP80870-2 100ug #AP80870-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMP-5,Bmp-5
Accession No.	Uniprot ID:P49003
Uniprot	P49003
GeneID	12160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bmp5

Note: This product is for in vitro research use only