

Bcl-2-binding component 3

Catalog No: #AP80872



Package Size: #AP80872-1 50ug #AP80872-2 100ug #AP80872-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bcl-2-binding component 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	JFY-1,p53 up-regulated modulator of apoptosis,PUMA
Accession No.	Uniprot ID:Q9BXH1
Uniprot	Q96PG8
GeneID	27113;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BBC3

Note: This product is for in vitro research use only