

Serine-protein kinase ATM

Catalog No: #AP80874



Package Size: #AP80874-1 50ug #AP80874-2 100ug #AP80874-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine-protein kinase ATM
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ataxia telangiectasia mutated,A-T mutated
Accession No.	Uniprot ID:Q13315
Uniprot	Q13315
GeneID	472;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATM

Note: This product is for in vitro research use only