

ATP synthase subunit beta, mitochondrial

Catalog No: #AP80875



Package Size: #AP80875-1 50ug #AP80875-2 100ug #AP80875-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ATP synthase subunit beta, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Accession No.	Uniprot ID:P99504
Uniprot	P99504
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP5B

Note: This product is for in vitro research use only