

# A disintegrin and metalloproteinase with thrombospondin motifs 7

Catalog No: #AP80876

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80876-1 50ug #AP80876-2 100ug #AP80876-3 1mg

## Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 7,COMPase
Accession No.	Uniprot ID:Q9UKP4
Uniprot	Q9UKP4
GeneID	11173;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADAMTS7

Note: This product is for in vitro research use only