

Axin-2

Catalog No: #AP80879



Package Size: #AP80879-1 50ug #AP80879-2 100ug #AP80879-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Axin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Axin-like protein,Axil,Axis inhibition protein 2,Conductin
Accession No.	Uniprot ID:O70240
Uniprot	O70240
GeneID	29134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Axin2

Note: This product is for in vitro research use only